



Technical Data Sheet

3-Ply Face Mask



Item # 91019P3

Contact Information:

Distributor: Dick Medical Supply, LLC.
Address: 664 N. Hague Avenue Suite B,
Columbus, Ohio 43204, USA
Phone: (614)-444-2300
Fax: (614)-444-2301
E-Mail: gkneeland@dms4ems.com
Revision: 05172020

664 N. Hague Avenue, Suite B, Columbus, Ohio 43204 (614) 444-2300 Office (614) 444-2301 Fax

Think DMS for EMS!



1. Product Overview

- Flat, easy on/easy off mask with standard pleats
- Easy to breathe through
- Bi-directional pleats for increased facial comfort and fit
- Has at least 95% of Bacterial Filtration Efficiency and Particulate Filtration Efficiency
- Meets American Society for testing and materials Level 1 standards

Flat, easy on/off ear loop mask with standard pleats. Lightweight, comfortable and easy to breathe through. Intended for use in environments with minimal fluid exposure. General patient care and other environments with minimum fluid exposure, isolation, IC practices, Central processing, SPD, clinics, long term care facilities and home health care

2. Specifications

Product Name:	3-Ply Earloop Procedure Face Mask
Color:	Blue
Size:	Universal
Standards:	Latex Free
Certification:	EN 149:2001 +A1:2009, FDA Registered, & CE Certified ASTM Level 1
Reusable:	No
Country of Origin:	China
Application:	Minimum fluid environments
Item#:	91019P3
Pack Size:	50ea/Box

Disclaimer: Dick Medical Supply, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. Dick Medical Supply makes no representations or warranties either express or implied, regarding the suitability of the material for any purpose or the accuracy of the information contained within this document. Accordingly, Dick Medical Supply will not be responsible for damages resulting from use of or reliance upon this information.