

High Output Nebulizer
which delivers medication
quickly and efficiently

A recent study compared the performance of 17 nebulizers regarding dead volume, time, MMAD, percent of aerosol mass in the respirable range, etc. The VixOne performed among the most efficient in each of the categories, and most importantly provided the largest respirable mass of aerosol available to the patient. — Chest, August 1996

Ideal particle size production (MMAD=2.7μ) allowing high lung deposition which results in maximum therapeutic effect.*

VixOne™ features and benefits:

- Anti-spill "Tee" design
- Leak-proof and secure seal with simple 1/4 turn assembly design
- Very low dead volume (0.7cc)
- High medication bowl capacity (10cc)
- Graduated volume (cc) markings on bowl
- Nebulizes in upright positions or at angles up to 45°
- Latex Free

* Reference CHEST Vol. 110, p.498-505 August 1996



Westmed
When Your Care is Critical



Product Number	Description	Qty
VixOne™ Nebulizer		
0214	VixOne, Nebulizer (Individually bagged)	50
VixOne™ Nebulizer with Accessories		
0200	VixOne, Multi Nebulizer Kit, Full Nebulizer Kit, 5 Extra Nebulizers	25
0209	VixOne, Nebulizer Mouthpiece, Universal "T", 7' Kink Resistant Tubing w/Threaded Nut & Flexhose	50
0210	VixOne, Nebulizer Mouthpiece, Universal "T", 7' Kink Resistant Tubing & Flexhose	50
0217	VixOne, Nebulizer Universal "T", 7' Kink Resistant Tubing (Unassembled)	50
0222	VixOne, Nebulizer Mouthpiece, Universal "T", 7' Kink Resistant Tubing	50
0248	VixOne, Nebulizer Pack 2 each VixOne Nebulizer (individually bagged), 2 each VixOne Nebulizer, Mouthpiece, Universal "T", 7' Kink Resistant Tubing, & Flexhose	25
0310	VixOne, Nebulizer w/7' Kink Resistant Tubing and Adult Mask	50
0311	VixOne, Nebulizer w/7' Kink Resistant Tubing and Pediatric Mask	50
0312	VixOne, Nebulizer w/7' Kink Resistant Tubing and Super Spike™ Pediatric Mask	50
0313	VixOne, Nebulizer w/7' Tubing, Threaded Nut and Super Spike™ Pediatric Mask	50
0319	VixOne, Nebulizer w/7' Tubing, Threaded Nut and Pediatric Mask	50
0326	Vix on Demand™, Patient Actuated Nebulizer	50
Infant Nebulizer Kits		
0235	VixOne Infant Extension Kit, w/ 15mm "T"	50
0236	VixOne Infant Nebulizer Kit, VixOne Nebulizer, Extension Set 3/8" Corrugated Tubing w/45° "T", 2" Flexible PVC Tube Adaptor, & 7' Kink Resistant Tubing	50
0238	Infant Extension Set, 3/8" Corrugated Tube with 45° "T", 2" Flexible PVC Tube Adaptor	50
VixOne™ Reusable Nebulizer		
0205	VixOne, Nebulizer (Individually bagged)	10
0254	VixOne, Nebulizer Mouthpiece, Universal "T", 7' Kink Resistant Tubing	10

VixOne™ is a Trademark of Westmed, Inc.

Aerosolized Medication Device for Infants. An alternative to aerosol "blow-by"

The PediNeb™ was designed to enhance the delivery of aerosolized medication, while reducing the apprehension encountered by young patients.

PediNeb™ features and benefits:

- Calms the infant/toddler during aerosol treatment. Medication is inhaled through the nose, not the mouth.
- Medication Port. Directs aerosolized medication to the nose and eliminates medication waste associated with the "blow-by" technique.
- VixOne small volume nebulizer. Clinically proven, providing superior medication delivery performance.
- Easy to operate, clean and efficient.
- Developed by a pediatric nurse practitioner.
- Latex free.

Distributed by:



The PediNeb™, a new modality in aerosolized medication delivery.

Patent Number 5,904,140
PediNeb™ is a Trademark of Westmed, Inc.

Product Number	Description	Qty
PediNeb™ Nebulizer		
0382	VixOne Nebulizer, Pacifier, 45 Degree Elbow, 12 & 7' Kink Resistant Tubing	
0383	VixOne Nebulizer, Pacifier, Pediatric Tee, 12 Poppel Extension Set & 7' Kink Resistant Tubing	
0385	Pedi-Neb, Pacifier with 45 Degree Elbow 25	