

C.H.A.R.L.I.E. Neonatal Resuscitation Simulator
product profile

Product Name	C.H.A.R.L.I.E. Neonatal Resuscitation Simulator
Addresses Training Needs/Curricula	<ul style="list-style-type: none">• Meets course curriculum for NRP® and Pediatric Advanced Life Support (PALS)• Suitable for pre- and in-hospital training (EMS professionals and nurses)
Configuration	<ul style="list-style-type: none">• Full-body manual simulator• Multiple Packages for Optimal Training• Heat-resistant (can be used in neonatal ICU unit and incubator)• Features two birth anomalies for simulation• SimVs offers an advanced simulation experience with C.H.A.R.L.I.E. NRP® Neonate with built-in scenario library and checklists
Dimension	<ul style="list-style-type: none">• 19 1/2" (49.53 cm) from crown to heel• Approximately 8 1/2 lbs (3.86 kg)
Product Code/SKU	LF01420 CHARLIE with Interactive ECG LF01421 CHARLIE LF0142101 CHARLIE with Sim VS Pre-Hospital Defibrillator Simulation System LF0142102 CHARLIE with Sim VS Transport Simulation System LF0142103 CHARLIE with Sim VS Nursing Essentials Simulation System LF0142104 CHARLIE with Sim VS Nursing Med-Surg Simulation System



Nasco Healthcare
16 Simulaids Drive
Saugerties, NY 12477
Phone: 1.833.NASCOHC (627-2642)
Email: info@nascohealthcare.com
www.nascohealthcare.com



C.H.A.R.L.I.E. Neonatal Resuscitation Simulator




Ready to save lives.
Right from the start.

C.H.A.R.L.I.E. Neonatal Resuscitation Simulator


Developed in conjunction with the latest Neonatal Resuscitation Program Guidelines (NRP®)

- C: Compression and cardiac
- H: Heat compatible can be put in an incubator without melting
- A: Airway, arterial access via the umbilicus
- R: Resuscitation
- L: Laryngeal mask capable for oxygen use or PPV
- I: Intravenous, intraosseous, intubation
- E: ECG



C.H.A.R.L.I.E. Neonatal Resuscitation Simulator:

Full-term newborn simulator for training on the most common neonatal care scenarios: Enables nurses and EMS workers to practice the essentials of neonatal care.



Add SimVs to upgrade to C.H.A.R.L.I.E. NRP® Neonate - an advanced simulation experience

- Includes an instructor's tablet that controls the displays of 3 types of monitor interfaces: defibrillator, hospital monitor, and mechanical defibrillator

This manikin is intended to be used as an instructional aide in the teaching and demonstrating of live rescue training. Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.