


SEMPERSHIELD® EXTENDED CUFF Nitrile EXAMINATION GLOVES

DISTRIBUTED BY
Life-Assist
Helping Heroes Save Lives



Extended Cuff Exam Nitrile Gloves

- **Thickness:** Competes with a 5-6 mil glove
- Extended cuff (min 280mm)
- Superior protection and strength
- Has been tested for protection against Fentanyl
- Excellent barrier and chemical resistance
- Meets USP 800 requirements for PPE compliance for chemotherapy
- Meets FDA 21 CFR Part 177 requirements for food handling

For more information about protection against fentanyl, please contact your sales rep.

100 / 10 /

REORDER INFO

Glove Color: Violet Blue

SKU	Size	Gloves/box	Boxes/case
SSEC102	Small	100 gloves/box	10 boxes/case
SSEC103	Medium	100 gloves/box	10 boxes/case
SSEC104	Large	100 gloves/box	10 boxes/case
SSEC105	X-Large	100 gloves/box	10 boxes/case

SAMPLES: Available by request - ask your sales representative

For maximum protection, we recommend changing nitrile gloves every 20-30 minutes

RECOMMENDED USES



PRODUCT SPECS



Rev07252022

sempermed®
YOUR GLOBAL GLOVE EXPERT™

100+
YEARS
- EST. 1920 -

Our Gloves are Disposable. The People they Protect are Not.®
www.sempermedusa.com • Ph 1.800.366.9545

